User Manual



FAVORGEN[®] FavorPrepTM Fungi/Yeast Genomic DNA Extraction Mini Kit

(For Research Use Only)

Kit Contents:

Cat. No:	FAFYG 000 (4 preps)	FAFYG 001 (50 preps)	FAFYG 001-1 (100 preps)
Beads Tube	4 pcs	50 pcs	100 pcs
FA Buffer	5 ml	60 ml	120 ml
FB Buffer	2.7 ml	32 ml	65 ml
TG1 Buffer	2 ml	27 ml	55 ml
TG2 Buffer	2 ml	15 ml	30 ml
W1 Buffer ^a (Concentrate)	1.3 ml	22 ml	44 ml
Wash Buffer ^b (Concentrate)	1 ml	10 ml	20 ml
Elution Buffer	0.5 ml	7 ml	15 ml
Lyticase Solution	250 µl	550 µl × 5	550 µl × 10
Proteinase K (Liquid)	100 µl	1050 µl	1050 µl × 2
TG Mini Column	4 pcs	50 pcs	50 pcs × 2
Collection Tube	8 pcs	100 pcs	100 pcs × 2
Elution Tube	4 pcs	50 pcs	100 pcs
User Manual	1	1	1
Preparation of W1 Buffer and Wash Buffer by adding ethanol (96~100%) and store at RT .			
Ethanol volume for W1 Buffer ^a	0.5 ml	8 ml	16 ml
Ethanol volume for Wash Buffer ^b	4 ml	40 ml	80 ml

Specification:

Principle: Sample size:	Mini spin column (silica matrix) 1~5×10° cells
Operation time:	<60 minutes
Column applicability:	Centrifugation and vaccum

Additional requirement to be provided by user

- 1. Microcentrifuge capable of speed at ~18,000 × g
- 2. 1.5 ml microcentrifuge tube
- 3.96~100% ethanol
- 4. Vortex mixer
- 5. Block heater or water bath

Important Notes:

- 1. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
- 2. Store the Lyticase Solution at -20°C on arrival.
- 3. Caution: Lyticase Solution and FB Buffer containing 14 mM of ß-mercaptoethanol is hazardous to human health. Perform the procedures involving Lyticase Solution and FB Buffer in a chemical fume hood.
- 4. Add required volume of ethanol (96~100%) to W1 Buffer and Wash Buffer for the first open. Store the buffer at room temperature.
- 5. Prepare block heater or water bath to 37°C for the step 4; 55°C for the step 9 before operation.

General Protocol:

Please Read Important Notes Before Starting Following Steps.

- 1. Transfer 1~5×10⁶ of cultured cells (fungi/yeast) to a 1.5 ml microcentrifuge tube. (not provided)
- 2. Add 1 ml of FA Buffer to the cells and resuspend the cells by pipetting.
- 3. Descend the cells by centrifuging at 5,000 × g for 2 minutes and discard the supernatant completely.
- 4. Resuspend the cells in 550 μl of FB buffer and add 50 μl of Lyticase Solution, mix well by vortexing. Incubate the sample at 37°C for 30 minutes.
 -Caution: Lyticase Solution and FB Buffer containing 14 mM of β-mercaptoethanol is hazardous to human health. perform the procedures involving Lyticase Solution and FB Buffer in a chemical fume hood.
- 5. (Optional) If RNA-free genomic DNA is required, add 8 µl of 50 mg/ml RNase A (not provided) and incubate for 2 minutes at room temperature.
- 6. Descend the cells by centrifuging at 5,000 × g for 10 minutes. Remove the supernatant completely.
- 7. Add 450 µl TG1 Buffer and mix well by pipetting. Transfer the sample mixture to a Bead Tube. (provided)
- Mix well by Plus-vortexing for 5 minutes.
 Elongate the pulse-vortexing time to 15~30 minutes if the sample cells are hard to be broken.
- 9. Add 20 µl of **Proteinase K and mix well by vortexing.** Incubate at 55°C for 15 minutes; vortex 30 seconds for every 5 minutes incubation.
- 10. Centrifuge the sample mixture at 5,000 × g for 1 minute and transfer 200 µl of the supernatant to a new 1.5 ml microcentrifuge tube (not provided).
- 11. Add 200 µl of TG2 Buffer and mix well by pipetting.
- 12. Add 200 µl of ethanol (96~100%) and mix well by pulse-vortexing for 10 seconds.
- 13. Place a **TG Mini Column** in **Collection Tube**. Transfer the sample mixture (including any precipitate) carefully to **TG Mini Column**. Centrifuge at 11,000 × g for 30 seconds **then place the TG Mini Column to a new Collection Tube**.
- Add 400 µl of W1 Buffer to the TG Mini Column. Centrifuge at 11,000 × g for 30 secconds. Discard the flow-through and then place the TG Mini Column back to the Collection Tube.
 -Make sure ethanol has been added into W1 Buffer at the first use.
- Add 750 µl of Wash Buffer to the TG Mini Column. Centrifuge at 11,000 × g for 30 secconds. Discard the flow-through and then place the TG Mini Column back to the Collection Tube.
 -Make sure ethanol has been added into Wash Buffer at the first use.
- Centrifuge at full speed (~18,000 × g) for an additional 3 minutes to dry the column.
 -Important Step! This step will avoid the residual liquid to inhibit subsequent enzymatic reactions.
- 17. Place the **TG Mini Column** to **Elution Tube**.
- Add 50~100 µl of Elution Buffer or ddH2O to the membrane center of the TG Mini Column. Stand TG Mini Column for 3 minutes.
 -Important Step! For effective elution, make sure that the elution solution is dispensed onto the membrane center and is absorbed completely.
- 19. Centrifuge at full speed (~18,000 × g) for 1 minute to elute total DNA.
- 20. Store total DNA at 4°C or -20°C.