

# FavorPrep™Total RNA Isolation Kit II

#### (For Research Use Only)

## **Kit Contents:**

Cat. No:	FATRS000 (5 Rxn)	FATRS050 (50 Rxn)	FATRS100 (100 Rxn)
Lysis Buffer	1.5 ml	12.5 ml	25 ml
2 M NaOAc, pH 5.2	150 µl	1.2 ml	2.5 ml
Wash Buffer (Concentrate) <sup>a</sup>	1 ml	7 ml	15 ml
Release Buffer	0.32 ml	3 ml	6 ml
RNA Column	5 pcs	50 pcs	100 pcs
Collection Tube	5 pcs	50 pcs	100 pcs
Preparation of Wash Buffer by adding	g ethanol (96~100%) and <b>s</b>	tore at RT.	
Ethanol volume for Wash Buffer <sup>a</sup>	4 ml	28 ml	60 ml

### **Specification:**

Principle:	mini spin column (silica matrix) <100 ug RNA/column
Sample size:	up to 1×10 <sup>6</sup> cultured cells
	up to 50 mg tissue
Operation time:	30 minutes
Column applicability:	centrifugation and vacuum

### Additional requirement to be provided by user

- 1. Microcentrifuge capable of speed at ~12,000 rpm
- 2. 1.5 ml microcentrifuge tube
- 3.96~100% ethanol
- 4. Water-Saturated Phenol
- 5. Chloroform
- 7. Vortex
- 8. Water bath or dry bath

### **Important Notes:**

- 1. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
- 2. Store the kit at room temperature.
- 3. Caution: phenol and chloroform are hazardous to human health. Perform the procedures involving phenol and chloroform in a chemical fume hood.
- 4. Add required volume of ethanol (96~100%) to Wash Buffer at the first open. Store the solution at room temperature.

### Protocol

Read the Important Note before starting the following steps.

Hint: Preheat Release Buffer to 65°C for step 12.

- 1. Add 200  $\mu$ l Lysis Buffer into the tube containing up to 100 mg tissue or 1×10<sup>6</sup> cultured cell pellet.
- 2. Vigorous mixing by vortexing. Incubate at room temperature for 10 minutes.
- 3. Add 20 µl 2 M NaOAc, pH 5.2.
- 4. Add 180  $\mu l$  water-saturated phenol and 40  $\mu l$  chloroform into the tube, vortex vigorously for 2 minutes.
- 5. Centrifuge at 12,000 rpm for 3 minutes. Transfer the upper phase into a clean tube.
- 6. Add 96~100% ethanol (2.33× volume of upper phase) to the upper phase and mix well by shaking vigorously.
- If the upper phase volume is 220 µl, add 513 µl of 96~100% ethanol to upper phase. The final ethanol concentration of whole mixture will be 70%.
- Note: Precipitates may be visible after addition of 96~100% ethanol. Resuspend precipitate completely by vigorous shaking and proceed immediately to step 7.
- 7. Transfer mixture, incoulding precipitate to RNA Column in the Collection Tube. Incubate for 1 minute.
- 8. Centrifuge at 12,000 rpm for 30 seconds.
- 9. Add 650  $\mu l$  Wash Buffer (ethanol added). Incubate for 1 minute.
- Make sure that ethanol has been added into Wash Buffer at the first open.
- 10. Centrifuge at 12,000 rpm for 1 minute to completely remove the residue liquid.
- 11. Put the RNA Column to a clean 1.5 ml tube.
- 12. Add 30~50 µl Release Buffer (preheated to 65°C) to the center of column. Incubate for 3 minutes.
- 13. Centrifuge at 12,000 rpm for 3 minutes to recover RNA.