



FavorPrep™ Soil DNA Isolation Mini Kit

Cat.No. : FASOI 000 (4 preps)
 FASOI 001 (50 Preps)
 FASOI 001-1 (100 Preps)

Kit Contents:

| | FASOI 000 (4 preps_sample) | FASOI 001 (50 preps) | FASOI 001-1 (100 preps) |
|-----------------------------|-------------------------------|-------------------------|----------------------------|
| Glass Beads | 4 vials | 50 vials | 100 vials |
| SDE1 Buffer | 3.6 ml | 40 ml | 70 ml |
| SDE2 Buffer | 1.2 ml | 15 ml | 25 ml |
| SDE3 Buffer | 1.2 ml | 15 ml | 30 ml |
| SDE4 Buffer | 1.5 ml | 25 ml | 40 ml |
| Wash Buffer (concentrate) * | 1.5 ml | 20 ml | 40 ml |
| Elution Buffer | 1.5 ml | 25ml | 50 ml |
| SDE Mini Column | 4 pcs | 50 pcs | 100 pcs |
| Collection Tube | 8 pcs | 100 pcs | 200 pcs |
| Elution Tube | 4 pcs | 50 pcs | 100 pcs |
| Bead tube | 4 pcs | 50 pcs | 100 pcs |
| User Manual | 1 | 1 | 1 |

* Preparation of Wash Buffer for first use:

| Cat. No: | FASOI 000 | FASOI 001 | FASOI 001-1 |
|--------------------------------|-----------|-----------|-------------|
| ethanol volume for Wash Buffer | 6 ml | 80 ml | 160 ml |

Description:

FavorPrep™ Soil DNA Isolation Mini Kit operates through our high-quality beads-beating disruption method and is perfect for use with diverse soil samples of up to 0.5g. Our silica membrane technology, and spin column along with beads-beating method guarantee the high-quality purification and isolation of DNA that can be used for PCR, genotyping, arrays, etc.

Specifications:

Principle: spin column (silica membrane)

Sample: 0.25 ~ 0.5 g

Operation time: < 60 min

Elution volume: 50 ~ 200 µl

Important Notes:


1. Buffers provided in this system contain irritants. Wear gloves, safety glasses, and a lab coat when handling these buffers.
2. Check SDE1 Buffer before use. Warm SDE1 Buffer at 60°C for 10 minutes if any precipitate forms.
3. Add indicated volume of ethanol (96-100%) to Wash Buffer before use.
4. Prepare a heating block or a water bath to 70 °C. If DNA is isolated from gram positive bacteria, prepare a heating block or a water bath to 95 °C for another incubation.
5. All centrifuge steps are done at full speed (~18,000 x g) in a microcentrifuge.
6. Preheat Elution Buffer or ddH₂O to 60°C for elution step.

General Protocol:

Please read Important Notes before starting with the following steps.

1. **Add the Glass Beads into a 2.0 ml Bead Tube (provided).**
 Transfer 0.25 ~ 0.5 g of soil sample into Bead Tube then place it on ice.
 --If the sample is liquid, add 200 µl of sample into a 2.0 ml Beads Tube.
2. **Add 600 µl of SDE1 Buffer to the sample, vortex at maximum speed for 5 minutes. Incubate the sample at 70 °C for 10 minutes and vortex the sample twice during the incubation.**
 --For isolation of DNA from gram positive bacteria, do a further incubation at 95 °C for 5 minutes.
3. **Briefly spin the tube to remove drops from the inside of the lid.**
4. **Cool down the sample mixture and add 200 µl of SDE2 Buffer. Mix well by vortexing. Incubate the sample on ice for 5 minutes.**
5. **Centrifuge at full speed (~ 18,000 x g) for 5 minutes.**
6. **Carefully transfer the clarified supernatant to a 1.5 ml microcentrifuge tube (not provided). Measure the volume of the supernatant.**
 --Avoid pipetting any debris and pellet.
7. **Add 1 volume of isopropanol and vortex to mix well. Centrifuge at full speed for 10 min to pellet DNA.**
 -- For example: If the clarified lysate volume is 450 µl, add 450 µl of isopropanol to the clarified lysate.
8. **Carefully discard the supernatant and invert the tube on the paper towel for 1 min to remove residual liquid.**
 --Do not disrupt the pellet.
9. **Add 200 µl of pre-heated Elution Buffer or ddH₂O, vortex to dissolve the DNA pellet completely.**
10. **Add 100 µl of SDE3 Buffer to the sample, mix well by vortexing. Incubate the sample at room temperature for 3 minutes.**
 --**Note:** SDE3 Buffer must be suspended completely by vigorously vortexing before every use.
 -- Cut off the end of 1 ml tip to make it easier to pipette the SDE3 Buffer.

Brief Procedure:



11. Centrifuge at full speed for 2 minutes.
12. Carefully transfer the supernatant to a 1.5 ml microcentrifuge (not provided), and measure the volume of the supernatant.
--Avoid pipetting any debris and pellet.
13. (Optional) If RNA-free DNA is required, add 1 µl of 100 mg/ml RNase A (not provided) to the sample and mix well. Incubate at room temperature for 2 min.
14. Briefly spin the tube to remove drops from the inside of the lid.
15. Add 1 volume of SDE4 Buffer and 1 volume of ethanol (96~100%). Mix thoroughly by pulse-vortexing.
For example: If the clarified lysate volume is 250 µl, add 250 µl of SDE4 Buffer and 250 µl of ethanol (96~100%) to the sample.
16. Place a SDE Column into a Collection Tube and transfer all of the sample mixture to the SDE Column. Centrifuge at full speed for 1 min, discard the flow-through, and place the SDE Column into a new Collection Tube.
17. Add 750 µl of Wash Buffer (ethanol added) to the SDE Column. Centrifuge at full speed for 1 min, then discard the flow-through.
--Make sure that ethanol (96~100%) has been added into the Wash Buffer upon first use.
18. Repeat step 17.
19. Centrifuge at full speed for an additional 3 min to dry the SDE column.
--Important step! This step will avoid the residual liquid to inhibit subsequent enzymatic reactions.
20. Place the SDE Column into a 1.5 ml microcentrifuge tube (not provided). Add 50 ~ 200 µl of preheated Elution Buffer or ddH₂O onto the membrane center of the SDE Column. Stand the SDE Column for 2 min at room temperature.
--Important step! For effective elution, make sure that the Elution Buffer or ddH₂O is dispensed onto the membrane center and is absorbed completely.
21. Centrifuge at full speed for 1 min to elute DNA.

Troubleshooting

| Problem | Possible reasons | Solutions |
|---|--|--|
| Low or no yield of genomic DNA | | |
| | Sample stored incorrectly | Store the stool sample at -20 °C. |
| | Low amount of cells in the sample | Increase the sample size |
| Poor cell lysis | | |
| | Poor cell lysis because of insufficient beads beating time | Extend the beads beating time. |
| Insufficient binding of DNA to column's membrane | | |
| | Ethanol is not added into sample lysate before DNA binding | Make sure that the correct volumes of ethanol (96- 100 %) are added into the sample lysate before DNA binding. |
| | Ethanol and sample lysate did not mix well before DNA binding | Make sure that Ethanol and sample lysate have been mixed completely before DNA binding |
| Incorrect preparation of Wash Buffer W1/W2 | | |
| | Ethanol is not added into Wash Buffer when first used | Make sure that the correct volumes of ethanol (96- 100 %) are added into the Wash Buffer upon first use. |
| | The volume or the percentage of ethanol is not correct for adding into the Wash Buffer | Make sure that the correct volumes of ethanol (96- 100 %) are added into the Wash Buffer upon first use. |
| Elution of DNA is not efficient | | |
| | pH of water (ddH ₂ O) for elution is acidic | Make sure the pH of ddH ₂ O is between 7.0-8.5. Use Elution Buffer (provided) for elution . |
| | Elution Buffer or ddH ₂ O is not completely absorbed by column membrane | After Elution Buffer or ddH ₂ O is added, stand the SD Column for 5 min before centrifugation. |
| Poor quality of genomic DNA | | |
| A260/A280 ratio of eluted DNA is low | Poor cell lysis | |
| | Poor cell lysis because of insufficient beads beating time | Extend the beads beating time. |
| A260/A280 ratio of eluted DNA is high | A lot of residual RNA in eluted DNA | Add 8 µl of RNase A (50 mg/ml) to the eluate and incubate at 37 °C for 10 minutes. After incubation, add 200 µl of SD2 Buffer and 200 µl of ethanol (96~100%), mix well by plus -vortexing. Then follow the general protocol starting from step 7. |