

# Material Safety Data Sheet

According to EC directive 2001/58/EC



## FavorPrep™ 96-well GEL/PCR Clean Up Purification Kit

**FAVORGEN**  
BIOTECH CORP.

Version 1

Processing date 12.23.2014

Printing date 12.23.2014

### 1. IDENTIFICATION OF THE SUBSTANCE/ PREPARATION

Commercial Product Name: FavorPrep™ 96-well GEL/PCR Clean Up Purification Kit

Cat. No.: FAPKE 000 (1 plate), FAPKE 001 (4 plates), FAPKE 002 (10 plates)

Supplier: Favorgen Biotech Corp.

Address: Ping-Tung Agricultural Biotechnology Park No.37, Nong-Ke Rd., Ping-Tung 908, Taiwan

Phone Number: 886-8-762-1829

Fax: 886-8-762-0791

Emergency telephone number: 886-8-762-1829

### 2. CHEMICAL CHARACTERISATION

Bind Buffer D1: 1 Piece (s)

Ingredients 25-50% M Guanidinium chloride; CAS-No.00050-01-1

Wash Buffer (concentrated): 1 Piece (s)

No hazardous substances in concentrations to be declared.

Elution Buffer: 1 Piece (s)

No hazardous substances in concentrations to be declared.

### 3. HAZARDS IDENTIFICATION

\*\*\*\*\* EMERGENCY OVERVIEW \*\*\*\*\*

Component 3, 4, and 5 – Occupational exposure presents little or no health hazard.

Component 1 and 2 - Warning! Irritant. Harmful if swallowed. May cause allergic skin reaction. Possible reproductive system hazard based on animal data.

\*\*\*\*\*

#### **Potential Health Effects:**

**Eye:** Can cause moderate irritation, tearing and reddening, but not likely to permanently injure eye tissue. Can cause severe irritation. Eye contact may result in corneal injury. Symptoms may include discomfort or pain, excess blinking and tear production, with marked redness and swelling of the conjunctiva. Temporary vision impairment (cloudy or blurred vision) is possible.

**Skin:** Can cause moderate skin irritation, defatting, and dermatitis. Not likely to cause permanent damage. May cause allergic skin reaction.

**Inhalation:** Can cause minor respiratory irritation, dizziness, weakness, fatigue, nausea, and headache. Can be corrosive to the respiratory tract causing severe irritation and tissue damage. If liquid mists are breathed into the lungs, they may be rapidly absorbed through the lungs and injure specific target organs. No toxicity expected from inhalation.

**Ingestion:** Mildly irritating to mouth, throat, and stomach. Can cause abdominal discomfort. Harmful if swallowed. May cause systemic poisoning. Corrosive to tissue. Can cause severe and permanent damage to mouth, throat, stomach. Aspiration may lead to lung damage. Ingestion of this product may result in central nervous system effects including

# Material Safety Data Sheet

According to EC directive 2001/58/EC



## FavorPrep™ 96-well GEL/PCR Clean Up Purification Kit

**FAVORGEN**  
BIOTECH CORP.

Version 1

Processing date 12.23.2014

Printing date 12.23.2014

headache, sleepiness, dizziness, slurred speech and blurred vision. Harmful if swallowed. May cause systemic poisoning.

Chronic: No data on cancer. Contains a substance(s) that is a possible reproductive system hazard based on high dose tests with laboratory animals.

### 4. FIRST AID MEASURES

**Eye:** Flush eyes with plenty of water for at least 20 minutes retracting eyelids often. Tilt the head to prevent chemical from transferring to the uncontaminated eye. Remove contact lenses, clean before re-use. Get immediate medical attention.

**Skin:** Wash with soap and water. Remove contaminated clothing and launder immediately, and discard contaminated leather goods, and wash before re-use. Get medical attention immediately if irritation develops or persists.

**Inhalation:** Remove to fresh air. If breathing is difficult, have a trained individual administer oxygen. If not breathing, give artificial respiration and have a trained individual administer oxygen. Get medical attention immediately.

**Ingestion:** **Complement 3, 4 and 5** – Give plenty of water, if conscious. If vomiting occurs naturally, wash mouth out. Be prepared to induce vomiting upon a physician's advice. Obtain medical attention if symptoms develop.

**Complement 1 and 2** - Do not induce vomiting and seek medical attention immediately. Drink two glasses of water or milk to dilute. Provide medical care provider with this MSDS. Corrosive.

Note To Physician: Treat symptomatically.

### 5. FIRE-FIGHTING MEASURES

Flashpoint Deg C: Not available.

Upper Flammable Limit %: Not available.

Lower Flammable Limit %: Not available.

Autoignition Temperature Deg C: Not available.

Extinguishing Media:

**Complement 3, 4 and 5** – Use means appropriate for surrounding materials.

**Complement 1 and 2** - Use alcohol resistant foam, carbon dioxide, dry chemical, or water spray when fighting fires. Water or foam may cause frothing if liquid is burning but it still may be a useful extinguishing agent if carefully applied to the fire. Do not direct a water stream directly into the hot burning liquid. Use water spray/fog for cooling.

**Firefighting Techniques/Equipment:** Do not enter fire area without proper protection including self-contained breathing apparatus and full protective equipment. Fight fire from a safe distance and a protected location due to the potential of hazardous vapors and decomposition products.

# Material Safety Data Sheet

According to EC directive 2001/58/EC



## FavorPrep™ 96-well GEL/PCR Clean Up Purification Kit

**FAVORGEN**  
BIOTECH CORP.

Version 1

Processing date 12.23.2014

Printing date 12.23.2014

Hazardous Combustion Products: Includes carbon dioxide, carbon monoxide, dense smoke.

### 6. ACCIDENTAL RELEASE MEASURES

Accidental releases may be subject to special reporting requirements and other regulatory mandates. Refer to Section 8 for personal protection equipment recommendations.

Spill Cleanup: **Complement 3, 4 and 5** – Absorb spill. Common absorbent materials should be effective.

Deposit in appropriate containers for removal and disposal.

**Complement 1 and 2** - Exposure to the spilled material may be irritating or harmful. Follow personal protective equipment recommendations found in Section VIII of this MSDS. Additional precautions may be necessary based on special circumstances created by the spill including; the material spilled, the quantity of the spill, the area in which the spill occurred. Also consider the expertise of employees in the area responding to the spill. Ventilate the contaminated area. Prevent the spread of any spill to minimize harm to human health and the environment if safe to do so. Wear complete and proper personal protective equipment following the recommendation of Section VIII at a minimum. Dike with suitable absorbent material like granulated clay. Gather and store in a sealed container pending a waste disposal evaluation.

### 7. HANDLING AND STORAGE

Storage of some materials is regulated by federal, state, and/or local laws.

Storage Pressure: Ambient

Handling Procedures:

**Complement 3, 4 and 5** – Keep closed or covered when not in use.

**Complement 1 and 2**– Harmful or irritating material. Avoid contacting and avoid breathing the material.

Use only in a well ventilated area. Keep closed or covered when not in use.

Storage Procedures: Suitable for most general chemical storage areas.

### 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Exposure Limits:

Component	OSHA PEL (ppm)	AGCIH TWA (ppm)
GUANIDINE HCL	Not established.	Not established.

Engineering Controls: Local exhaust ventilation or other engineering controls are normally required when handling or using this product to avoid overexposure.

Personal Protective Equipment:

Eye: An eye wash station must be available where this product is used. Wear chemical goggles.

Skin: Avoid skin contact by wearing chemically resistant gloves, an apron and other protective equipment depending upon conditions of use. Inspect gloves for chemical break-through and replace at regular intervals. Clean protective equipment regularly. Wash hands and other exposed areas

# Material Safety Data Sheet

According to EC directive 2001/58/EC



## FavorPrep™ 96-well GEL/PCR Clean Up Purification Kit

**FAVORGEN**  
BIOTECH CORP.

Version 1

Processing date 12.23.2014

Printing date 12.23.2014

with mild soap and water before eating, drinking, and when leaving work. Have a safety shower available.

Respiratory: A respiratory protection program that meets OSHA's 29 CFR 1910.134 and ANSI Z88.2 requirements must be followed whenever workplace conditions warrant a respirator's use.

### 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance/physical state: Liquid solution / suspension

Order: No odor.

Specific Gravity/Density: Not established.

Octanol/water Partition Coeff: Not established.

Volatiles: Not established.

Evaporation Rate: Not established.

Viscosity: Not established.

### 10. STABILITY AND REACTIVITY

Stability: Stable under normal conditions.

Conditions to Avoid: Strong oxidizing agents. High temperatures. Strong acids.

Hazardous Decomposition Products: None known.

Hazardous Polymerization: Hazardous polymerization will not occur.

### 11. TOXICOLOGICAL INFORMATION

Acute Toxicity:

Inhalation/Respiratory: GUANIDINE: 475 MG/KG 3250 MG/KG

Target Organs: Kidneys. Central Nervous System.

Carcinogenicity:

NTP: Not tested.

IARC: Not listed.

OSHA: Not regulated.

Other Toxicological Information

Contains a substance(s) that is a possible reproductive system hazard based on high dose tests with laboratory animals.

Toxicity (others): LD50 intraperitoneal mouse 500 mg/kg

### 12. ECOLOGICAL INFORMATION

Ecotoxicological Information: No ecological information available.

Environmental Fate (Degradation, Transformation, and Persistence): Bioconcentration is not expected to occur.

### 13. DISPOSAL CONSIDERATIONS

# Material Safety Data Sheet

According to EC directive 2001/58/EC



## FavorPrep™ 96-well GEL/PCR Clean Up Purification Kit

**FAVORGEN**<sup>®</sup>  
BIOTECH CORP.

Version 1

Processing date 12.23.2014

Printing date 12.23.2014

Regulatory Information: Not applicable.

Disposal Method: Clean up and dispose of waste in accordance with all federal, state, and or local environmental.

### 14. TRANSPORT INFORMATION

Proper Shipping Name: Not Determined.

Subsidiary Hazards:

### 15. OTHER INFORMATION

WHMIS Rating 0-4:

FIRE:0.

HEALTH: 2.

REACTIVITY: 0

#### Abbreviations

N/A - Data is not applicable or not available

SARA - Superfund and Reauthorization Act

HMIS - Hazard Material Information System

WHMIS - Workplace Hazard Materials Information System

NTP - National Toxicology Program

OSHA - Occupational Health and Safety Administration

IARC - International Agency for Research on Cancer

PROP 65 - California Safe Drinking Water and Toxic Enforcement Act of 1986

EINECS - European Inventory of Existing Commercial

Chemical Substances