


Nucleic Acids Molecular weight markers



20bp DNA Ladder, ready-to-use (20bp- 500bp)

50bp DNA Ladder (50bp -500pb)

100bp DNA Ladder (100bp – 3,000bp)


250bp DNA Ladder (250bp – 5,000bp)

1Kb DNA Ladder (500bp-10,000bp)


1Kb DNA Ladder Plus (500bp -10,000bp)

Lambda Hind III Marker / Lambda Hind III Plus Marker (0.125 -23.130Kb)


Favorgen® DNA Markers




20bp DNA Ladder, ready-to-use







50bp DNA Ladder



100bp DNA Ladder, ready-to-use



250bp DNA Ladder



1Kb DNA Ladder

1Kb DNA Ladder Plus

Lambda Hind III Marker

Lambda Hind III Plus Marker

Order Information

Cat. No.	Product Name	Conc.	Size Range	Bands	Loading	Q' ty	Package
FALA20	Favorgen® 20bp DNA Ladder	400 ng/1μl	20bp-500bp	25	2.5μl	100	250μl

FALA50	Favorgen® 50bp DNA Ladder	50 ng/1µl	50bp-500bp	10	2.5 – 5.0µl	100-125	500µl
FALA100	Favorgen® 100bp DNA Ladder ready-to-use	140 ng/1µl	100-3000bp	13	5.0µl	100	500µl
FALA101	Favorgen® 100bp DNA Ladder	100 ng/1µl	100bp-3000bp	13	5.0µl	100	500µl
FALA250	Favorgen® 250bp DNA Ladder	200 ng/1µl	250bp-5,000bp	13	2.5µl	100	250µl
FALA200	Favorgen® 1Kb DNA Ladder	250 ng/1µl	500bp-10,000bp	10	2.0µl	200	400µl
FALA201	Favorgen® 1Kb DNA Ladder Plus	250 ng/1µl	250bp-10,000bp	12	2.0µl	200	400µl
FALH100	Favorgen® Lambda Hind III marker	250ng/1µl	0.13-23.130Kb	8	2.5µl	200	500µl
FALH101	Favorgen® 1Kb DNA Ladder Plus	150 ng/1µl	0.13-23.130Kb	10	2.5µl	200	500µl