


Protein Standards


FavorStain™ Prestained Markers

Prestained Protein Markers for SDS-PAGE and Western Blotting provide a quick and easy way to identify and characterize the separation of molecules during electrophoresis and to assess transfer efficiency on blots.



Favorgen® Unstained Protein Markers

Favorgen Unstained Protein Markers allow accurate molecular weight determination on SDS-PAGE gels. They are available for SDS-PAGE, IEF, 2-D and Western Blotting.



Order Information

Cat. No.	Product Name	Concen.	Size Range	Bands	Q' ty	Package
FAPM-002	Favorgen® Unstained Protein Marker, low range	0.2-0.3mg/ml	14, 21, 31, 45, 66, 97kDA	6	100	250µl x 2 vials
FAPM-002-1	Favorgen® Unstained Protein Marker, low range	0.2-0.5mg/ml	3.4, 6.5, 14, 21, 31, 45, 66, 97kDA	8	100	250µl x 2 vials
FAPM-001	FavorPrep Prestained Protein Marker, low range (single color)	3.0mg/ml	14, 24, 31, 47, 76, 100kDA	8	100	250µl x 2 vials
FAPM-001-1	FavorPrep Prestained Protein Marker, low range (dual color)	2.6mg/ml	14, 24, 31, 54, 76, 97kDA	8	100	250µl x 2 vials
FAPM-001-2	FavorPrep Prestained Protein Marker, low range (multicolor)	6.8mg/ml	14, 21, 31, 54, 76, 97kDA	8	100	250µl x 2 vials